

2702 Three Channel Cordless Falls Monitor

cura¹
Helping To Care



Monitor up to 3
Cordless Devices



2702



Nurse Call

... cordless with economy

2702 Three Channel Cordless Falls Monitor

Combining the most essential features of the 2700 but with greater economy and simplicity while still maintaining professional standards. Ideal for home care or caring for a single patient.

- Monitor up to 3 cordless devices simultaneously.
- Comes with batteries.
- Monitor auto re-arm if the patient gets back into bed.
- Amber LED visual alarm indicator.
- Battery operated or optional external plug pack (2406).
- Up to 6 months battery life (nominal).
- Easy to access battery compartment.
- Can be connected to a nurse call system.
- 3 selectable volume levels.
- Alerts for lost nurse call connection.
- Fault supervision for pressure pads (low battery, device lost).



WHAT'S IN THE BOX?



COMPATIBLE PRODUCTS

- All cura¹ cordless sensors & buttons.
- 3504 Mounting Bracket
- 2742 Soft Monitor Cover
- 2406 Plug Pack 6 Vdc



cura¹
Helping To Care

10 Gray Street, Kilkenny, South Australia 5009
sales@cura1.com | www.cura1.com
+61 (0) 8 8268 8343 (International)

1300 12-CARE (2273)