

Three Channel Cordless Falls Monitor



PART NUMBER: 2702

FEATURES

- Modern ergonomic design
- Monitor up to 3 cordless devices simultaneously
- NO / NC switch in battery compartment
- Comes with batteries
- Alarm auto re-arm
- Amber LED visual alarm indicator
- Battery operated or optional external plug pack (2406)
- Up to 6 months battery life (nominal)
- Can be connected to a nurse call system
- 3 selectable volume levels
- Alerts for lost nurse call connection, lost pad, lost AC power supply & low battery
- Alerts for low battery in paired bed / chair pad transmitters

TECHNICAL SPECIFICATIONS

- DC operation using either a 3 x AA batteries (included) or 6 V plug pack (2406)
- Dimensions 120 x 83 x 35 mm
- Weight 210 g without batteries
- Up to 50 m* range when used with cordless bed pads and chair pads
- Operating temperature 5-40° C, relative humidity < 85%
- Operating frequency 433 Mhz fully supervised
- Voltage free normally open / normally closed relay contact alarm output
- Low battery detection
- Housing: ABS (UL-94HB)

ACCESSORIES

All Cordless Bed & Chair Pads & Floor Mats (transmitters 2719, 2720, 2732, 2670), 2721 iBeam, 3503 Wireless Door Contact, 2595 Wristband Alarm Transmitter, 2599 Wireless SOS Call Button, 3502 Wireless Handset, 2558 Wireless Call button with Reset & Pull Cord, 2588 PET, 2698 Universal Transmitter, 2742 Soft Monitor Cover, 3504 ActiveCare Mounting Bracket, 2643 Nurse Call Cable, 2406 Plug Pack 6 Vdc

cura¹.com

*The range of wireless devices may vary and be dependent on your specific environmental conditions.
cura¹ products are designed to augment well designed care programs and practices, cura¹ does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.