

HERO+

Hero Medical Wheelchair Range

*Choose Your Hero: Comfortable,
Functional and Stylish Wheelchairs for a
Wide Range of Needs and Budgets*



Selecting the Right Wheelchair For You

The right wheelchair can provide independence and comfort for those experiencing mobility challenges. Finding the best wheelchair to meet a user's needs is essential to maximising the benefit that a wheelchair can provide.

There are a few key considerations when selecting a wheelchair:



Type

There are two main types of wheelchairs, each catering to different needs: transit and self-propelled wheelchairs.

Transit wheelchairs, also known as attendant propelled wheelchairs, are generally lighter in weight, have smaller rear wheels and designed to be pushed from behind.

Self-Propelled wheelchairs give users the option to push themselves or be pushed by an attendant.

Comfort

It is important to consider the comfort of both the user and the attendant when selecting a wheelchair.

A variety of features can contribute to comfort, including a padded seat and backrest, cushioned armrests, and height-adjustable leg rests that can easily be removed as needed.

Also, push handles with ergonomic hand grips aid the attendant in pushing the chair with ease.



Style

Some users place an importance on the appearance of a wheelchair. The Hero Medical range of wheelchairs has options that are both functional and stylish.

Weight

Different settings may have requirements that impact wheelchair weight. Community or facility settings may require sturdier frames, which are often heavier.

Lightweight frames are often easier to transport, such as lifting in and out of a car.

Premium Comfy Lite

The Hero Medical Comfy Lite Wheelchair Range is a premium, lightweight wheelchair range with ergonomic design and enhanced safety features.

Compact Ideal for storing and transporting with removable and folding parts

Lightweight Perfect for portability and handling due to its lightweight aluminium frame

Ergonomic Splayed arms, soft grooved arm pads and padded seat and back covers provide maximum comfort

Stylish Matte black finish with removable, two-tone seat and back covers

Safety* Premium functions include lockable attendant brakes, side release seat belt and telescopic anti-tippers



Type
Transit (AP)



Self-Propelled (SP)

Colour ● Matte Black

Seat Width 16" / 18" / 19"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S309152	AP	Matte Black	40cm (16")	48cm	93cm	45cm	13"	9.7kg	135kg
S309153	AP	Matte Black	45cm (18")	48cm	93cm	45cm	13"	9.9kg	135kg
S309154	AP	Matte Black	47.5cm (19")	48cm	93cm	45cm	13"	10.1kg	135kg
S309156	SP	Matte Black	40cm (16")	48cm	93cm	45cm	22"	11.7kg	135kg
S309157	SP	Matte Black	45cm (18")	48cm	93cm	45cm	22"	11.9kg	135kg
S309158	SP	Matte Black	47.5cm (19")	48cm	93cm	45cm	22"	12.1kg	135kg

*Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing. Compliant to ISO 7176-19 certification - Crash Testing (Please note: S309152 & S309153 are crash test approved).

Active Ultra Lite

The Hero Medical Active Ultra Lite Wheelchair is a premium, stylish, ultra-lightweight wheelchair, designed for easy transportation and optimum user comfort.

Compact Ideal for storing and transporting with removable and folding parts

Lightweight Perfect for portability and handling due to it's lightweight aluminium frame


Ergonomic Height adjustable, curved armrests and dual layer cushioned seat and backrest promote user comfort

Stylish Finished with a metallic carbon texture finish for a sophisticated look

Safety* Premium functions, including lockable attendant brakes, side release seat belt and anti-tippers



Type  Transit (AP)  Self-Propelled (SP)

Colour  Carbon Texture

Seat Width 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S309155	AP	Carbon Texture	46cm (18")	50cm	93cm	44cm	12"	8.9kg	135kg
S338462	SP	Carbon Texture	46cm (18")	50cm	93cm	44cm	22"	11.2kg	135kg

*Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing.

Feather Lite

The Hero Medical Feather Lite Wheelchair is an ultra-light, stylish, and transport friendly wheelchair that is great value for money.

- Compact** Ideal for storing and transporting with removable and folding parts
- Lightweight** Perfect for portability and handling due to it's lightweight aluminium frame
- Great Value** An affordable option with many features included
- Safety*** Features lockable attendant brakes, seat belt and anti-tippers



Type  Transit (AP) **Colour**  Blue  Red  Silver **Seat Width** 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S301725	AP	Blue	46cm (18")	47.5cm	91cm	40.6cm	8"	8.4kg	115kg
S301726	AP	Red	46cm (18")	47.5cm	91cm	40.6cm	8"	8.4kg	115kg
S301724	AP	Silver	46cm (18")	47.5cm	91cm	40.6cm	8"	8.4kg	115kg

*Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing.

Super Budget

The Hero Medical Super Budget Wheelchair is a sturdy, entry level wheelchair that is great value for money.

- Sturdy** Steel frame is sturdy, robust, and highly durable
- Easy to Clean** Vinyl material enables effective hygiene and infection control - ideal for facility settings
- Great Value** An affordable option with many features included
- Safety*** Features park brakes and seat belt



Type  Transit (AP)  Self-Propelled (SP) **Colour**  Glossy Black **Seat Width** 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S310065	AP	Glossy Black	45cm (18")	48cm	92cm	42.5cm	13"	15.4kg	115kg
S310064	SP	Glossy Black	45cm (18")	48cm	92cm	42.5cm	24"	17.6kg	115kg

*Compliant to AS/NZS 3695.1:2011 Standards – Strength, Durability, Load, Static, Impact and Fatigue Testing.



Feature Comparison

	Super Budget		Feather Lite	Active Ultra Lite		Comfy Lite	
	<i>Transit</i>	<i>Self-Propelled</i>	<i>Transit</i>	<i>Transit</i>	<i>Self-Propelled</i>	<i>Transit</i>	<i>Self-Propelled</i>
Structural Details							
Frame Colour	Glossy Black	Glossy Black	Blue, Red, Silver	Carbon Texture	Carbon Texture	Matte Black	Matte Black
Frame Material	Steel	Steel	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
Collapsible Frame	+	+	+	+	+	+	+
Seat							
Removable Cover						+	+
Material	Vinyl	Vinyl	Fabric	Fabric	Fabric	Fabric	Fabric
Seat Belt	+	+	+	+	+	+	+
Wheels							
Quick Release					+		+
Solid Tyres	+	+	+	+	+	+	+
Rear Wheel Size	13"	24"	8"	12"	22"	13"	22"
Front Castor Wheel Size	8"	8"	6"	8"	8"	6"	6"
Anti Tipper Wheels			+	+	+	+	+
Legrests							
Height Adjustable	+	+	+	+	+	+	+
Foot Strap	+	+	+	+	+	+	+
Swing-away	+	+	+	+	+	+	+
Removable	+	+	+	+	+	+	+
Armrests							
Height Adjustable				+	+		
Padded	+	+	+	+	+	+	+
Backrests							
Half Folding Backrest			+	+	+	+	+
Removable Cover						+	+
Backrest Pocket	+	+	+	+	+	+	+
Material	Vinyl	Vinyl	Fabric	Fabric	Fabric	Fabric	Fabric
Brakes							
Attendant Brakes			+	+	+	+	+
User Accessible Brakes	+	+			+	+	+
Other Features							
Standards Testing (ASNZ or ISO)	+	+	+	+	+	+	+
Safe Work Load (SWL)	115kg	115kg	115kg	135kg	135kg	135kg	135kg
Net Weight	15.4 kg	17.6 kg	8.9kg	8.9kg	11.2kg	9.7-10.1kg	11.7-12.1kg