

ROYALE

ROBOTOER

Precision mobility for life on the move

● E60-A

● E60 Pro-A



The Royale Robooter E60-A and E60 Pro-A Powerchairs are built for versatility – compact enough for indoor use, capable enough for outdoors. Their omnidirectional wheels and foldable frame make them the perfect companion for daily routines, providing stability, comfort and confidence across varied terrain.

E60-A E60 Pro-A

All-terrain powerchairs balancing agility and comfort

Designed with ergonomic comfort and modern convenience in mind, they're ready for everyday adventures.

The E60-A and E60 Pro-A make outdoor independence easy. Their omnidirectional wheels and shock-absorbing suspension handle pathways and parks with confidence, while the portable, foldable frame ensures convenient transport and storage.

E60-A
Grey Frame
10" wheels



E60 Pro-A
Green Frame
12" wheels



Foldable frame

Compact for transport



Reclining backrest

Customised comfort at 8°–38°



Omnidirectional wheels

Exceptional manoeuvrability



Complementary headrest included



E60-A vs **E60 Pro-A**

- Grey frame
- 10" omnidirectional wheels
- Suitable for indoor and outdoor use
- Green frame
- 12" omnidirectional wheels
- Suitable for more rugged terrain


Attendant control

Via app or bluetooth remote



Storage included

Rear backpack & underseat basket



Coil-spring suspension

Smooth, comfortable ride



Specifications

E60-A **E60 Pro-A**

Frame Colour	Grey	Green
Frame Material	Aluminium	Aluminium
Safe Working Load (SWL)	150 kg	150 kg
Battery Type	Lithium-Ion	Lithium-Ion
Battery Size	Standard 20 AH Optional: 12 AH or 25 AH	Standard: 20 AH Optional: 12 AH or 25 AH
Range	17 km (12 AH) 23 km (20 AH) 29 km (25 AH)	17 km (12 AH) 23 km (20 AH) 29 km (25 AH)
Max Speed	8 km/h	10 km/h
Size (W x L x H Assembled)	64 x 98.7 x 93.2 cm	64 x 105.2 x 95.6 cm
Seat Width	42 cm	42 cm

E60-A **E60 Pro-A**

Seat Depth	43 cm	43 cm
Backrest Height	46 cm	46 cm
Backrest Angle	8° – 38°	8° – 38°
Distance Between Armrests	45.5 – 51.5 cm	45.5 – 51.5 cm
Weight (Without Battery)	30.8 kg	34 kg
Front Wheels	10" omni	12" omni
Rear Wheels	10" pneumatic	12" pneumatic
Max Safe Incline Slope	9°	15°
Turning Circle	80 cm	85 cm
Clearance	6 cm	8 cm
Folded Size (W x L x H)	64 x 77 x 50 cm	64 x 85.5 x 52.5 cm