

ROYALE

ROBOTOOTER

Effortless mobility for everyday freedom

E80



The Royale Robooter E80 Powerchair brings together carbon fibre engineering and contemporary design for mobility on the go. Weighing just 15.9 kg, it's one of the lightest powered wheelchairs in its class – built for convenience, style and comfort.

E80

Stylish carbon fibre powerchair built for life on the go

The ideal choice for those who want the freedom to go anywhere – without the weight.

Every detail of the E80 is made for effortless mobility. The carbon fibre frame balances strength and style while the intuitive folding feature makes transport and storage easy. Flexible attendant controls and LED lighting add smart functionality.



Carbon fibre frame

Just 15.9kg for easy lifting



Foldable frame

Compact for transport



Ergonomic seat angle

Supports posture





E80

- ▶ **Sleek gradient finish**
Modern aesthetic
- ▶ **LED front & rear lights**
Improved visibility
- ▶ **Ergonomic seat & backrest**
Comfort for daily use

Attendant control


Via app or bluetooth remote

Storage included

Rear backpack & underseat basket

Airline-compatible battery

12 AH battery is ideal for travel



Specifications

Frame Colour	Black / White
Frame Material	Carbon Fibre
Safe Working Load (SWL)	150 kg
Battery Type	Lithium-Ion
Battery Size	Standard: 12 AH Optional: 20 AH or 25 AH
Range	18 km (12 AH) 25 km (20 AH) 29 km (25 AH)
Max Speed	8 km/h
Size (W x L x H Assembled)	62.8 x 99.8 x 96.5 cm
Seat Width	42 cm

Seat Depth	43 cm
Backrest Height	51 cm
Backrest Angle	18°
Distance Between Armrests	48 cm
Weight (Without Battery)	15.9 kg
Front Wheels	8" PU
Rear Wheels	10" pneumatic
Max Safe Incline Slope	6°
Turning Circle	85 cm
Clearance	4 cm
Folded Size (W x L x H)	62.8 x 71.5 x 44.5 cm