



GET AMONGST IT!



QUANTUM® OUTBACK® 4x4
shown in Tan with 10" (25.4cm)
safe seat elevation and 14"
(38.1cm) pneumatic tires



POWER BEYOND THE PAVEMENT

The Quantum® Outback® 4x4 is the ultimate power chair for off-road independence, built to tackle rugged terrain with confidence. Equipped with powerful 4-pole sealed motors, all-wheel drive, and 14" pneumatic knobby tyres, it delivers outstanding traction and control. Advanced SRS (Smooth Ride Suspension) and robust construction ensure a stable, comfortable ride over gravel paths, bush tracks, and uneven ground. Whether you're exploring country properties or outdoor trails, the Outback® 4x4 offers the freedom to go further combining power, durability, and smart design for life beyond the pavement.



QUANTUM®



quantumrehab.com.au

Specifications

DRIVE WHEELS	38.1cm (14") knobby Pneumatic Tyres
DRIVETRAIN	4 high torque motors
SUSPENSION	Automotive-grade smooth ride suspension
MAXIMUM SPEED ¹	Up to 10 km/h (5.6 km/h when elevated)
MAXIMUM OBSTACLE CLIMBING ABILITY	10.2 cm
MAXIMUM RATED SLOPE	8° (14.1%)
GROUND CLEARANCE	10.7 cm (at motors)
TURNING RADIUS ²	125 cm
OVERALL LENGTH ²	106.7 cm
BASE WIDTH	73.2 cm
SEATING SIZES	
TRU-Balance® 3 Static	Width: 16-24"(40.6cm - 61cm) Depth: 16-24"(40.6cm - 61cm)
TRU-Balance® 3 Power Tilt	Width: 16-24"(40.6cm - 61cm) Depth: 16-24"(40.6cm - 61cm)
TRU-Balance® 3 Power Lift & Tilt	Width: 16-24"(40.6cm - 61cm) Depth: 16-22"(40.6cm - 55.9cm)
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Power Tilt ^{5,9}	19-21.5" (48.3 cm - 54.5 cm)
TRU-Balance® 3 Power Lift & Tilt ^{5,9}	20-21" (50.8 cm - 53.4 cm)
BATTERY SIZE	
	MK 88Ah(2)
BATTERY WEIGHT ⁷	
	28 kg each
RANGE PER CHARGE ³	
	Up to 25.75 km at 136 kg Up to 27.4 km at 99.8 kg
AVAILABLE ELECTRONICS	
	120A Q-Logic 3
BATTERY CHARGER	
	8A Off-board
MOTOR PACKAGES	
	4-pole motors
WEIGHT CAPACITY ⁴	
	160 kg
BASE WEIGHT	
	226 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 1.5kg. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded.

Standard Features

- 4-pole sealed motors
- Maximum speed of up to 10km/h 5.6 km/h when elevated)
- LED lighting package
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions

Options

- 10" (25.4 cm) of power adjustable safe seat elevation at speeds of up to 5.6 km/h
- Quantum® backup camera
- Glove box
- Backpack holder
- Multi-use accessory bag
- Hydration system
- Accessory bar

Outdoor utility mount

Colors



Tan (Matte)

Sage (Matte)

Camo Black (Matte)

Camo White (Matte)

Camo dessert (Matte)

