

# Q6 edge<sup>®</sup> 2.0



## ENGINEERED FOR CONFIDENT PERFORMANCE

*The Q6 Edge<sup>®</sup> 2.0 serves as Quantum's high-performance mid-wheel drive chair ideal as a dependable daily mobility companion. With standard high-speed (10 km/h) motors and ATX suspension, it delivers responsive control and stability—especially effective when custom spring packages are specified for lighter or paediatric clients.*

*Optional iLevel<sup>®</sup> seat elevation (up to 30.5 cm / 12") at walking speed further enhances daily independence. The Edge 2.0 supports a wide range of seating and electronics options, offering exceptional value and versatility.*



| [quantumrehab.com.au](http://quantumrehab.com.au)

## Specifications

<b>DRIVE WHEELS (Mid-Wheel Drive)</b>	35.5 cm or (14") Pneumatic or Flat Free Tyres 14x4" wide Flat Free Tyres
<b>CASTER WHEELS</b>	<b>Front:</b> 15.2 cm (6") solid <b>Rear:</b> 15.2 cm (6") solid
<b>ANTI -TRIP WHEELS</b>	N/A
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 10 km/h Up to 5.6 km/h w/ iLevel activated
<b>GROUND CLEARANCE</b>	6.8 cm (2.7") battery tray
<b>TURNING RADIUS</b>	52 cm (20.5") w/ 15.2 cm (6") castors 59 cm (23.2") w/ 20.3 cm (8") castors
<b>OVERALL LENGTH <sup>2</sup></b>	90 cm (35.5") w/ 15.2 cm (6") castors 102.6 cm (40.4") w/ 20.3 cm (8") castors
<b>BASE WIDTH</b>	61.9 cm (24.4")
<b>SEATING SIZES ( TRU-BALANCE® 3 )</b>	
<b>TRU-Balance 4®</b>	N/A
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 30.5cm - 55.8cm (12-22") <b>Depth:</b> 30.5cm - 55.8cm (12-22")
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 30.5cm - 55.8cm (12-22") <b>Depth:</b> 30.5cm - 55.8cm (12-22")
<b>TRU-Balance® 3 Power Lift &amp; Tilt</b>	<b>Width:</b> 30.5cm - 55.8cm (12-22") <b>Depth:</b> 30.5cm - 55.8cm (12-22")
<b>SEAT-TO-FLOOR HEIGHTS ( TRU-BALANCE® 3 )</b>	
<b>TRU-Balance® 4</b>	N/A
<b>TRU-Balance® 3 Static Seat</b>	42cm - 47cm (16.5-18.5")
<b>TRU-Balance® 3 Power Tilt <sup>5/7</sup></b>	44.5cm - 49.5cm (17.5-19.5")
<b>TRU-Balance® 3 Tilt &amp; iLevel® <sup>5/7</sup></b>	44.5cm - 49.5cm (17.5-19.5")
<b>BATTERY SIZE</b>	MK 50Ah(2)
<b>BATTERY WEIGHT <sup>6</sup></b>	17.2 kg each
<b>AVAILABLE ELECTRONICS</b>	120A Q-Logic 3
<b>RANGE PER CHARGE <sup>3</sup></b>	Up to 32 km
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-pole
<b>WEIGHT CAPACITY <sup>4</sup></b>	136 kg
<b>BASE WEIGHT</b>	193 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.  
Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 1.5kg. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded.

## Standard Features

- 2-pole motors
- ATX Suspension provides stability for smooth obstacle transitions
- Q-Logic 3 Drive controls
- Speeds up to 10 kph
- LED front and rear fender lights
- USB charger included
- ISO 7176-19 transit crash tested approved with all configurations

## Options

- TRU-Balance® 3 Power Positioning:  
iLevel® option allows 30.5 cm 12" of seat elevation while at walking speed (5.6 km/h)
- Power Tilt, Recline, and Elevating Legrests that can power to Floor for easy transfers.
- Weight specific spring packages available
- Accu -Trac advanced tracking technology
- Optional 8" castors all round
- NEW backup camera
- LED Lighting Package

## Colors

