

PRIDE®

MOBILITY

Take your
outdoor adventures
to the next level.



Introducing our new hard top canopy.

Whether cruising downtown, exploring new trails, or just enjoying a day outdoors, protect yourself from the elements with our new hard top canopy.

Compatible with the Celebrity MFP® Outback®, Pathrider® 130XL - 140XL and Pathrider® Endurance, its durable roof and windshield are perfect for the sun, wind, and other weather you may encounter throughout your journeys. Take your outdoor adventures to the next level with this exciting accessory.



Available on these selected models



PATHRIDER 130 XL
PATHRIDER 140 XL

PATHRIDER®
— ENDURANCE —

OPTIONAL canopy weather cover (not shown) with integrated vents will protect your product and canopy when not in use.

OPTIONAL jumbo rear basket

OPTIONAL tool clips (pair)

Shown on Outback



Rugged ABS plastic roof for superb durability and impact resistance

OPTIONAL Large panoramic rearview mirror

Strong polycarbonate windshield

Durable, powder-coated steel frame

For precise alignment and efficient application of the canopy, a two-person installation is recommended.

Optional extra accessories.



Tool clips (pair)
(OPTIONAL ACCESSORY)



Jumbo rear basket
(OPTIONAL ACCESSORY)



Canopy weather cover
(OPTIONAL ACCESSORY)

PRIDE[®]
MOBILITY