

ROYALE

Over Chair & Bed Swivel Table

Dual-purpose support and tabletop for everyday living



The Royale Over Chair Swivel Table combines practical functionality with a sleek, modern design. Dual swivel movement and a height-adjustable support handle provide stable support, flexible positioning and a convenient surface within easy reach.

Over Chair & Bed Swivel Table

Easy height adjustment for confident support

Simple hydraulic push-button allows smooth, tool-free height adjustment, making it easy to position the handle at a comfortable and supportive height.

Dual swivel mechanism for flexible positioning

Support arm allows the tabletop to swing smoothly over the individual for use, then pivot away to provide clear access when standing or transferring.

Adjustable base for safer compatibility

Adjustable base designed to slide easily under a wide range of chairs, recliners and lounges. The base length can be extended to suit different furniture depths, helping achieve a stable, secure fit beneath the seating.

Stylish tabletop with everyday practicality


Laminated tabletop available in two contemporary finishes – Black Ebony or Brown Mahogany – and contains a built-in cup holder and utensil slot to keep drinks and small essentials close to hand.



Ergonomic support handle with button to adjust height



Tabletop with 360° arm swivel & 360° axis swivel



Adjustable base width & non-slip rubber base pads

Tabletop Colour Options



Black Ebony
S339491



Brown Mahogany
S339492

Specifications

Maximum support handle load	130 kg
Maximum tabletop load	9 kg
Tabletop pivot	360°
Tabletop rotation	360°
Tabletop size	60 x 40 cm
Height adjustment (tabletop)	70 – 104.5 cm
Height adjustment (support handle)	85 – 119.5 cm
Base feet adjustment (inner length)	49 – 81 cm