

# Presenting the Cobalt Rotating Bed

**COBALT**  
HEALTH

## Advanced support

The Cobalt Rotating Bed is designed to promote independence and provide advanced levels of support, comfort and mobility for individuals who need some assistance getting in and out of bed.

The profiling chair function sees the entire bed surface rotate smoothly, providing full-body support and minimising both user and carer strain, whilst reaching a sitting position for easier bed transfers.



### Cobalt Rotating Bed

Code: B26-RB

#### PRODUCT FEATURES

Package includes Rotating bed, head & footboards, fold down side rails, mattress and mattress cover (sheet)

Electronic height adjustment range 35cm – 70cm (not including mattress)

Functions include profiling chair, height, backrest, knee bend, Trendelenburg and Anti-Trendelenburg

Optional 20cm bed length extension, including bolster and cover



## Increased independence

The Rotating Bed can assist users to maintain their independence by supporting them to get in and out of bed without the need for any additional family or caregiver intervention.

Where caregiver assistance is required the Rotating Bed can help to reduce strain on the carer by electronically reaching the sitting position and making transfers easier and safer.

Bed dimensions 98 x 215cm  
Mattress dimensions 90 x 190 x 8cm

Homestyle wooden beech head and footboards

Optional fabric slip on head & footboard covers in blue or grey

SWL 185kg (max user weight 165kg)


Warranty Frame 5 years, Electrics 2 Years, Mattress 1 Year



Fabric bed end covers



Replacement cover sheet



Foot end extension



Extension bolster

HELPING  
CARE  
HAPPEN

FOR FURTHER INFORMATION OR DEMONSTRATIONS, PLEASE CONTACT US

#### Victoria

[vicsales@cobalthhealth.com.au](mailto:vicsales@cobalthhealth.com.au)

(03) 8521 8444

43 Paramount Bvd, Cranbourne West VIC 3977  
1800 417 822

#### New South Wales

[sales@cobalthhealth.com.au](mailto:sales@cobalthhealth.com.au)

(02) 8853 1100

34/9 Salisbury Rd, Castle Hil NSW2154  
Cobalthhealth.com.au