

Manual wheelchairs

Panthera Micro 3



permobil

Smooth moves for little adventurers

The Panthera Micro 3 is designed for very young users and intended for indoor environments and supervised outdoor activities on flat surfaces. This lightweight, easy-to-handle manual wheelchair enables young children to explore the world around them.

Key benefits:

- Easy to push and manoeuvre
- Robust and lightweight
- Accessories for individual comfort



Specifications

	Micro 3	Micro 3 Long
Frame material	Chrome molybdenum tubing	Chrome molybdenum tubing
Seat width	24 - 27cm	24 - 27cm
Seat depth	15 - 20cm	20 - 25cm
Seat height	Front and rear 30.5cm (18" rear wheels) Front and rear 33cm (20" rear wheels)	Front and rear 30.5cm (18" rear wheels) Front and rear 33cm (20" rear wheels)
Seat angle	0°	0°
Total length	63cm	71cm
Total width	46 - 49cm	46 - 49cm
Backrest height	22 - 28cm	22 - 28cm
Backrest angle	90°	90°
Footplate distance	9.5 - 24.5cm	9.5 - 24.5cm
Camber	10°	10°
Transport weight	2.2 - 2.4kg	2.3 - 2.4kg
Wheelchair weight	3.4 - 3.5kg	3.5 - 3.6kg
Maximum user weight	30kg	30kg

Personalise your ride*

ANODIC BLACK

TELE-MAGENTA

RUBY RED

PEARL VIOLET

PINE GREEN

Permobil reserves the right to change specifications and design without prior notice. *More colour options available.

Actual product colours may vary.