

HERO+

Hero Medical Outdoor Lite All-Terrain Seat Walker

Take on outdoor terrain with confidence,
control and a comfortable place to rest.



Hero Medical Outdoor Lite All-Terrain Seat Walker

The Hero Medical Outdoor Lite All-Terrain Seat Walker is built for confident movement across outdoor environments. Large pneumatic wheels and rear suspension support a smoother, more controlled ride over uneven surfaces, while the lightweight frame remains easy to handle and transport. With a padded seat and adjustable backrest, it also provides a comfortable place to rest, supporting longer outings with greater ease.

Features

- **All-terrain pneumatic wheels:** support smooth movement across grass, gravel and uneven paths
- **Rear-wheel suspension:** reduces impact for improved comfort outdoors
- **Wide, stable frame:** enhances balance and confidence on varied terrain
- **Padded seat with adjustable backrest:** provides a comfortable resting option when needed
- **Height-adjustable handles:** support upright posture and a personalised fit
- **Light-touch handbrakes with parking lock:** provide controlled stopping and added safety



Pneumatic wheels with suspension



Removable storage bag

Specifications

Safe Working Load	136 kg
Frame Weight	8.5 kg
Gross Weight	9.9 kg
Seat Height	56 cm
Seat Width	46 cm
Handle Height	86 - 98.5 cm
Outer Width	64 cm
Front Castor	12.5"
Rear Castor	10"
Frame Material	Aluminium