

HERO+

Hero Medical Premium Knee Walker

Move confidently through recovery with
smooth, stable mobility support



Hero Medical Premium Knee Walker

The Hero Medical Premium Knee Walker offers a practical, more comfortable alternative to crutches during lower limb recovery. A padded, adjustable knee platform supports positioning, while rear-wheel suspension helps deliver a smoother, more controlled ride across indoor and outdoor surfaces. Its lightweight, compact frame makes it easy to manoeuvre in tight spaces, supporting day-to-day mobility while recovery progresses.

Features

- **Padded depth-adjustable knee platform:** supports comfort and positioning during extended use
- **Height-adjustable design:** adjustable handlebars and knee platform allows a personalised fit
- **Rear-wheel suspension:** reduces impact for a smoother, more controlled ride
- **Stable four-wheel base:** enhances balance and confidence when moving
- **Lightweight, compact frame:** easier to manoeuvre, transport and store
- **Responsive hand brakes:** provide controlled stopping for added safety



Depth-adjustable knee pad



Rear-wheel suspension

Specifications

Safe Working Load	136 kg
Overall Height	81-105 cm
Overall Length	86 cm
Overall Width	41 cm
Knee Pad Depth	34.5 - 38.5 cm
Knee Pad Width	17.5 cm
Knee Pad Height	43 - 59 cm
Handle Height	81 - 105 cm
Front Castor	9.5"
Rear Castor	8"