

# HERO+

## Hero Medical Wheelchair Range

*Choose Your Hero: Comfortable,  
Functional and Stylish Wheelchairs for a  
Wide Range of Needs and Budgets*



# Selecting the Right Wheelchair For You

The right wheelchair can provide independence and comfort for those experiencing mobility challenges. Finding the best wheelchair to meet a user's needs is essential to maximising the benefit that a wheelchair can provide.

There are a few key considerations when selecting a wheelchair:



## Type

There are two main types of wheelchairs, each catering to different needs: transit and self-propelled wheelchairs.

Transit wheelchairs, also known as attendant propelled wheelchairs, are generally lighter in weight, have smaller rear wheels and designed to be pushed from behind.

Self-Propelled wheelchairs give users the option to push themselves or be pushed by an attendant.

## Comfort

It is important to consider the comfort of both the user and the attendant when selecting a wheelchair.

A variety of features can contribute to comfort, including a padded seat and backrest, cushioned armrests, and height-adjustable leg rests that can easily be removed as needed.

Also, push handles with ergonomic hand grips aid the attendant in pushing the chair with ease.



## Style

Some users place an importance on the appearance of a wheelchair. The Hero Medical range of wheelchairs has options that are both functional and stylish.

## Weight

Different settings may have requirements that impact wheelchair weight. Community or facility settings may require sturdier frames, which are often heavier.

Lightweight frames are often easier to transport, such as lifting in and out of a car.

# Carbon Lite

The Hero Medical Carbon Lite is a premium lightweight wheelchair that sets the standard. Carbon fibre construction achieves remarkable weight reduction while delivering standout aesthetics and superior comfort.

**Compact** Folds flat for effortless transport and storage

**Carbon Fibre** The lightest frame in the range. Carbon fibre construction delivers exceptional weight reduction for easier lifting and handling

**Comfortable** Padded seat, backrest and calf pad provide superior comfort for longer outings

**Stylish** Modern carbon fibre aesthetic with a premium, contemporary finish

**Practical** Removable components make transport, storage and cleaning effortless

**Safety\*** Attendant brakes, anti-tip wheels and seat belt for confident, safer mobility



**Type**  Transit (AP)  Self-Propelled (SP)

**Colour**  Red & Black

**Seat Width** 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S339481	AP	Red & Black	48 cm (18")	48 cm	95 cm	50 cm	10.6"	8.8 kg	130 kg
S339480	SP	Red & Black	48 cm (18")	48 cm	95 cm	50 cm	22"	9.8 kg	130 kg

\*Compliant to AS/NZS 3965.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing.

# Premium Comfy Lite

The Hero Medical Comfy Lite is a feature-rich wheelchair range designed for everyday use – available in multiple sizes to suit a broader range of users, with superior ergonomic comfort and enhanced safety for longer sitting periods.

**Compact** Folds compactly with removable parts for straightforward storage and transport

**Lightweight** Aluminium frame keeps the chair manageable for everyday lifting and handling

**Ergonomic** Splayed arms, grooved arm pads and generous seat and back padding deliver maximum comfort for extended sitting periods

**Stylish** Sleek matte black finish with removable two-tone covers for a polished, practical look

**Safety\*** Lockable attendant brakes, side-release seat belt and telescopic anti-tippers – plus crash-tested models available for added peace of mind



**Type**  Transit (AP)  Self-Propelled (SP) **Colour**  Matte Black **Seat Width** 16" / 18" / 19"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S309152	AP	Matte Black	40 cm (16")	48 cm	93 cm	45 cm	13"	9.7 kg	135 kg
S309153	AP	Matte Black	45 cm (18")	48 cm	93 cm	45 cm	13"	9.9 kg	135 kg
S309154	AP	Matte Black	47.5 cm (19")	48 cm	93 cm	45 cm	13"	10.1 kg	135 kg
S309156	SP	Matte Black	40cm (16")	48 cm	93 cm	45 cm	22"	11.7 kg	135 kg
S309157	SP	Matte Black	45 cm (18")	48 cm	93 cm	45 cm	22"	11.9 kg	135 kg
S309158	SP	Matte Black	47.5 cm (19")	48 cm	93 cm	45 cm	22"	12.1 kg	135 kg

\*Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing. Compliant to ISO 7176-19 certification - Crash Testing (Please note: S309152 & S309153 are crash test approved).

# Active Ultra Lite

The Hero Medical Active Ultra Lite is a lightweight, versatile wheelchair designed for easy transportation and comfortable use across a range of everyday settings.

**Compact** Removable and folding parts make it quick to pack, store and transport

**Lightweight** Lightweight aluminium frame makes lifting in and out of vehicles easy - ideal for users on the go

**Ergonomic** Height-adjustable curved armrests and dual-layer cushioned seat and backrest adapt to a range of users and settings

**Stylish** Distinctive metallic carbon texture finish for a sophisticated look

**Safety\*** Lockable attendant brakes, side-release seat belt and anti-tippers for reliable everyday safety



**Type**  Transit (AP)  Self-Propelled (SP)

**Colour**  Carbon Texture

**Seat Width** 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S309155	AP	Carbon Texture	46cm (18")	50 cm	93 cm	44 cm	12"	8.9 kg	135kg
S338462	SP	Carbon Texture	46cm (18")	50 cm	93 cm	44 cm	22"	11.2kg	135kg

\*Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing.

# Feather Lite

The Hero Medical Feather Lite Wheelchair is an ultra-light, stylish, and transport friendly wheelchair that is great value for money.

**Compact** Ideal for storing and transporting with removable and folding parts

**Lightweight** Perfect for portability and handling due to it's lightweight aluminium frame

**Great Value** An affordable option with many features included

**Safety\*** Features lockable attendant brakes, seat belt and anti-tippers



Type  Transit (AP)

Colour  Blue  Red  Silver

Seat Width 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S301725	AP	Blue	46cm (18")	47.5cm	91cm	40.6cm	8"	8.4kg	115kg
S301726	AP	Red	46cm (18")	47.5cm	91cm	40.6cm	8"	8.4kg	115kg
S301724	AP	Silver	46cm (18")	47.5cm	91cm	40.6cm	8"	8.4kg	115kg

\*Compliant to AS/NZS 3695.1:2011 certification - Strength, Durability, Load, Static, Impact and Fatigue Testing.

# Super Budget

The Hero Medical Super Budget Wheelchair is a sturdy, entry level wheelchair that is great value for money.

- Sturdy** Steel frame is sturdy, robust, and highly durable
- Easy to Clean** Vinyl material enables effective hygiene and infection control - ideal for facility settings
- Great Value** An affordable option with many features included
- Safety\*** Features park brakes and seat belt



**Type**  Transit (AP)  Self-Propelled (SP) **Colour**  Glossy Black **Seat Width** 18"

Code	Type	Colour	Seat Width	Seat Height	Overall Height	Backrest Height	Rear Wheels	Net Weight	SWL
S310065	AP	Glossy Black	45cm (18")	48cm	92cm	42.5cm	13"	15.4kg	115kg
S310064	SP	Glossy Black	45cm (18")	48cm	92cm	42.5cm	24"	17.6kg	115kg

\*Compliant to AS/NZS 3695.1:2011 Standards – Strength, Durability, Load, Static, Impact and Fatigue Testing.



## Feature Comparison

	Super Budget		Feather Lite	Active Ultra Lite		Comfy Lite		Carbon Lite	
	Transit	Self-Propelled	Transit	Transit	Self-Propelled	Transit	Self-Propelled	Transit	Self-Propelled
<b>Structural Details</b>									
Frame Colour	Glossy Black	Glossy Black	Blue, Red or Silver	Carbon Texture	Carbon Texture	Matte Black	Matte Black	Red & Black	Red & Black
Frame Material	Steel	Steel	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Carbon Fibre	Carbon Fibre
Collapsible Frame	+	+	+	+	+	+	+	+	+
<b>Seat</b>									
Removable Cover						+	+	+	+
Material	Vinyl	Vinyl	Fabric	Fabric	Fabric	Fabric	Fabric	Fabric	Fabric
Seat Belt	+	+	+	+	+	+	+	+	+
<b>Wheels</b>									
Quick Release					+		+		+
Solid Tyres	+	+	+	+	+	+	+	+	+
Rear Wheel Size	13"	24"	8"	12"	22"	13"	22"	10.6"	22"
Front Castor Wheel Size	8"	8"	6"	8"	8"	6"	6"	6.5"	6.5"
Anti Tipper Wheels			+	+	+	+	+	+	+
<b>Legrests</b>									
Height Adjustable	+	+	+	+	+	+	+	+	+
Foot Strap	+	+	+	+	+	+	+	+	+
Swing-away	+	+	+	+	+	+	+	+	+
Removable	+	+	+	+	+	+	+	+	+
<b>Armrests</b>									
Height Adjustable				+	+				
Padded	+	+	+	+	+	+	+	+	+
<b>Backrests</b>									
Half Folding Backrest			+	+	+	+	+	+	+
Removable Cover						+	+	+	+
Backrest Pocket	+	+	+	+	+	+	+	+	+
Material	Vinyl	Vinyl	Fabric	Fabric	Fabric	Fabric	Fabric	Fabric	Fabric
<b>Brakes</b>									
Attendant Brakes			+	+	+	+	+	+	+
User Accessible Brakes	+	+			+	+	+		+
<b>Other Features</b>									
Standards Testing (ASNZ or ISO)	+	+	+	+	+	+	+	+	+
Safe Work Load (SWL)	115kg	115kg	115kg	135kg	135kg	135kg	135kg	130kg	130kg
Net Weight	15.4kg	17.6kg	8.9kg	8.9kg	11.2kg	9.7-10.1kg	11.7-12.1kg	8.8kg	9.8kg