

Finding the Right Robooter

ROBOTEER	E40	E80	E60-A	E60 Pro-A
Frame & Wheels				
Frame Colour	White	White/Black Gradient	Grey	Green
Frame Material	Aluminium	Carbon Fibre	Aluminium	Aluminium
Weight (without battery)	23 kg	15.9 kg	30.8 kg	34 kg
Front Wheels	8" Solid	8" Solid	10" Omni	12" Omni
Rear Wheels	10" Pneumatic	10" Pneumatic	10" Pneumatic	12" Pneumatic
Turning Circle	90 cm	85 cm	80 cm	85 cm
Seating				
Backrest Recline	Fixed 18°	Fixed 18°	8° - 38°	8° - 38°
Performance				
Max Speed	8 km/h	8 km/h	8 km/h	10 km/h
Range*	17-29 km	18-29 km	17-29 km	17-29 km
Max Incline	6°	6°	9°	15°
Safe Working Load	150 kg	150 kg	150 kg	150 kg
Features				
Airline Suitability*	✓ (with smaller 12AH battery option installed)	✓	✗	✗
Wireless Attendant Control	✓	✓	✓	✓

*Range & airline suitability varies depending on battery size (12 AH, 20 AH or 25 AH). See individual product specifications on our website for details.

ROYALE ROBOTEER Powerchair Range

The Royale Robooter range redefines what a powerchair can be – stylish, intelligent and built for real life.

From compact indoor manoeuvrability to confident outdoor performance, these folding powerchairs bring together smart technology, modern design and everyday comfort.

ILS INDEPENDENT LIVING SPECIALISTS

Explore the Royale Robooter range and book a free trial today

ilsau.com.au | 1300 008 267



PD1074-0326



With four models to choose from, there's a Robooter for every lifestyle – whether you're travelling, moving through daily routines or exploring the outdoors.

E40



Ideal for
Home & Community Use

Lightweight and compact. Perfect for apartment living and short community trips.

- 23 kg aluminium frame for easy carer handling
- 90 cm turning circle for tight spaces
- Flip-up armrests for safe transfers

E80



Ideal for
Active Travellers

ULTRA-Light carbon fibre frame built for easy day-to-day travel & airline friendly portability.

- 15.9 kg - lightest in the range
- 85 cm turning circle
- Airline-compatible battery for travel

E60-A



Ideal for
Everyday Versatility

Versatile performance for both indoor and outdoor environments.

- 80 cm turning circle for tight indoor spaces
- Omnidirectional wheels and suspension for agility and comfort
- Adjustable recline (8°–38°) for postural comfort

E60 Pro-A



Ideal for
Outdoor & All-Terrain Use

All-Terrain capability with enhanced suspension and 12" wheels for stability and control outdoors.

- 12" wheels for all-terrain traction
- 85 cm turning circle with large omnidirectional wheels
- 15° incline rating for outdoor confidence

Built for comfort and confidence, every Robooter features intuitive joystick control, a wireless attendant control remote included, bluetooth connectivity and integrated LED lighting for added safety.

Foldable frames make transport and storage easy, while thoughtful ergonomic details support all-day comfort.