

ROYALE

Air Carbon Fibre Folding Scooter

Ultralight carbon fibre construction and travel-ready design for confident use in the community, on holidays and everywhere in between.



Confident mobility, made portable

The Royale Air Carbon Fibre Folding Scooter delivers reliable powered mobility in an ultralight, foldable design that's built for life on the move.

Ultralight at just 13.7 kg

The carbon fibre frame keeps overall weight low, making it simple to lift, transport and store.

Tows like luggage when folded

A quick-fold mechanism and compact footprint make it easy to move through airports, cruise terminals and busy travel spaces.

Two batteries, up to 16 km range

Dual interchangeable lithium batteries keep you going longer, whether it's a local errand or a full day out.



Airline-friendly batteries*

Dual 6 AH batteries are ideal for travel

Compact footprint

Easy to store in car boots or tight spaces

Personalised comfort

Adjustable tiller and removable cushion for comfort and support



*Prior to travelling, please contact your airline to ensure your scooter and its battery comply with current airline policy.



- ▶ **Ultralight carbon fibre**
Modern finish and just 13.7 kg (without battery)
- ▶ **Quick-fold mechanism**
Folds in seconds and tows like luggage
- ▶ **NFC card start system**
Keyless start with controlled access and added security

Flip-up armrests

Assist safer transfers on and off the scooter



Front and rear LED lights

For improved visibility in low light conditions



Clear LCD display

Easy-to-read driving data at a glance



Specifications

Safe Working Load	120 kg	Assembled Size (W x L x H)	51.7 x 103 x 92 cm
Battery Type	Lithium	Folded Size (W x L x H)	52.5 x 44.5 x 69 cm
Battery Size	6 AH (x2)	Weight	13.7 kg (without battery) Each battery weighs 0.8 kg
Range	16 km (2 batteries) 8 km (1 battery)	Seat Width	47 cm
Maximum Speed	6 km/h	Seat Depth	38 cm
Maximum Safe Incline	8°	Seat Height	53 cm
Turning Circle	130 cm	Display	Digital
Ground Clearance	2 cm	Front Wheels	7"
Obstacle Climbing	5 cm	Rear Wheels	7"