

HERO+

Carbon Lite Wheelchair

Designed for easier lifting, compact storage and everyday mobility support, this range brings a modern, lightweight approach to manual wheelchairs.



Designed for lighter everyday mobility

A modern manual wheelchair range combining lightweight carbon fibre construction with compact folding designs for easier transport, storage and daily mobility support.

Compact folding frames, half-folding backrests and lightweight construction make the Carbon Lite range easier to lift into a vehicle boot, store at home or take away when travelling. Quick-release rear wheels on the Self-Propelled model further reduce transport weight.

Removable swing-away footrests, padded seating, anti-tip wheels and practical safety features support comfortable everyday mobility.



Carbon fibre frame

Lightweight construction with a sleek, contemporary finish



Compact folding design

Folds down easily for transport and storage



Removable rear wheels

Self-propelled model features push-pin release rear wheels



Two models for different mobility needs

Transit

Compact rear wheels for easier attendant manoeuvrability



Self-propelled

Supports independent mobility with 22" rear wheels



Removable footrests

Supports safer transfers and easier positioning



Padded seating

Everyday comfort for regular mobility use



Practical safety features

Anti-tip wheels, brakes and seat belt included



Specifications

Type	Self-Propelled	Transit
Frame Material	Carbon Fibre	Carbon Fibre
Colour	Red & Black	Red & Black
Safe Working Load (SWL)	130 kg	130 kg
Net Weight	9.8 kg	8.8 kg
Frame Weight	6 kg	6 kg
Seat Width	48 cm	48 cm
Seat Depth	42 cm	42 cm
Seat Height (floor to seat)	48 cm	48 cm
Overall Height	95 cm	95 cm
Backrest Height (from seat)	50 cm	50 cm
Armrest Height (from seat)	22 cm	22 cm
Front Castor	6.5"	6.5"
Rear Wheel	22"	10.6"